
































1

**Assay Data Sheet - Expected Ranges
BioProfile® FLEX Controls - LEVEL 1**

1

Parameter	Units	Expected Range		Exp. Date: 2018-01
		Lower Limit	Upper Limit	
Lot Number:	 619446			
pH	-	 6.641	 6.741	 6.616
PO ₂	mm Hg	 58.2	 70.6	 55.1
PCO ₂	mm Hg	 54.3	 67.9	 50.9
Gluc	g/L	 0.75	 1.25	 0.62
Lac	g/L	 0.78	 1.28	 0.65
NH ₄ ⁺	mmol/L	 0.61	 1.11	 0.48
Na ⁺	mmol/L	 69.7	 85.3	
K ⁺	mmol/L	 1.71	 2.29	
Ca ⁺⁺	mmol/L	 0.50	 0.74	
PO ₄	mmol/L	 3.41	 4.61	
p	mg/L	 105.7	 142.9	
Osmolality	mOsm/kg H ₂ O	 218	 228	



200 Prospect Street • Waltham, MA U.S.A. 02454-9141 • (781) 864-0800

LFN 08239B 2011-03