


































**3**

**Assay Data Sheet - Expected Ranges  
BioProfile® FLEX Controls - LEVEL 3**

**3**

Lot Number:  514239		Exp. Date: 2016-11	
		+/- 3 SD Expected Range	
Parameter	Units	Lower Limit	Upper Limit
pH	-	 7.318	 7.418
PO <sub>2</sub>	mm Hg	 134.2	 162.8
PCO <sub>2</sub>	mm Hg	 13.8	 20.2
Gluc	g/L	 12.75	 17.25
Lac	g/L	 4.25	 5.75
NH <sub>4</sub> <sup>+</sup>	mmol/L	 8.24	 11.24
Na <sup>+</sup>	mmol/L	 166.7	 203.3
K <sup>+</sup>	mmol/L	 5.55	 6.85
Ca <sup>++</sup>	mmol/L	 2.25	 2.85
Osmolality	mOsm/kg H <sub>2</sub> O	 603	 619
		Lower Limit	Upper Limit
		 7.293	 7.443
		 127.1	 169.9
		 12.2	 21.8
		 11.62	 18.38
		 3.87	 6.13
		 7.49	 11.99

*nova*  
biomedical

200 Prospect Street • Waltham, MA U.S.A. 02454-9141 • (781) 864-0800

LFN 08241A 2007-08